



Replacement Saline Concentrate

Mix with clean, potable water to create a preserved, buffered, saline solution. For use with the Fendall Porta Stream I, II, III



Overview

Key Features

Mix with clean, potable water (10:1 ratio) to create preserved, buffered and physiologically correct saline solution. Does not require costly plumbing. Provides a constant fluid flow rate and stream height for the entire 15-minute flushing period required by ANSI Z358.1 when used with Porta Stream II and III.

Hazards

- Airborne Particulates
- Biohazard
- Chemical
- Contamination
- Gas, Vapors, Smoke
- Splash

Regulations

- ANSI Z358.1 - ANSI Z358.1

Historical Brand

Fend-All

Warranty Information

LIMITED WARRANTY. Except for any warranty for a particular product or part specifically stated in our product literature, we warrant all goods sold under this Agreement against defects in materials and workmanship for one year from shipment. In the event of such defects, we will, at our sole option and as your exclusive remedy, repair or replace the defective product without charge, or refund the price you have paid. In order for this warranty to be effective, you must notify us within 15 days of discovery of the defect. You must return the defective product to us within the warranty period, properly packaged and with insurance and transportation costs prepaid. We will reship the product to you at our expense in the United States. This warranty does not cover defects which we determine are because of accident, alteration, modification, negligence, misuse, abuse, failure to perform adequate maintenance, normal wear and tear and normal maintenance. WE MAKE NO OTHER

Replacement Saline Concentrate

WARRANTIES, EXPRESS OR IMPLIED, AND WE SPECIFICALLY EXCLUDE ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NO PERSON IS AUTHORIZED TO EXTEND ADDITIONAL WARRANTIES ON OUR BEHALF.

Specifications

Solution	Saline Concentrate - Mix with clean, potable water (10:1 ratio) to create preserved, buffered and physiologically correct saline solution.
Expiration	30 month shelf life when unopened / 6 months once opened, from date of manufacture
Station Refilling	Prior to refilling the station must be cleaned to remove contamination from the station walls.

SKUs & Ordering Information

SKUs	Details
32-000513-0000	Sperian Saline Concentrate 180 oz. Saline Concentrate for Fendall Porta Stream I, II or other 14-16 gallon stations Case Qty: 4 Case Dimentions: 19" x 14" x 10" Case Weight (lbs): 52
32-000513-0020	Sperian Saline Concentrate - English/French (Canada) 180 oz Saline Concentrate Eng/French for Fendall Porta Stream II, III or Other 14-16 Gallon Stations Case Qty: 4 Case Dimentions: 19"x14"x10" Case Weight (lbs): 52
32-000509-0000	Sperian Saline Concentrate 70 oz. Saline Concentrate for Fendall Porta Stream I Case Qty: 6 Case Dimentions: 13" x 11" x 10" Case Weight (lbs): 32
32-000509-0020	Sperian Saline Concentrate - English/French (Canada) 70 oz. Saline Concentrate for Fendall Porta Stream I (English/French) Case Qty: 6 Case Dimentions: 13"x11"x10" Case Weight (lbs): 32

Replacement Saline Concentrate

Training & Proper Use

- Manuals provided for installation.
- Must be mixed with clean potable water.
- Prior to refilling the station must be cleaned to remove contamination from the station walls.

Technical Support

US Sperian Customer Care: 866-257-7374

RELATED PRODUCTS



For general inquiries
800.343.3411

[View other Sperian Sites](#)
© 2010 Sperian Protection